

NEW IN THE RANGE

SPIDER VEINS & VASCULAR · NEW

VasoClear 980 980 NM VASCULAR LASER

Spider veins gone. Skin left calm.

980 nm diode laser for vascular lesions

980

nm WAVELENGTH

30 W

LASER POWER

3

SPOT SIZES



KEY BENEFITS

- 980 nm is absorbed by oxyhaemoglobin: it targets the vessel with minimal melanin absorption
- Facial spider veins (nose, cheeks, chin), broken capillaries and leg thread veins
- Up to 90% clearance of fine spider veins in one session (manufacturer data)
- No anaesthetic cream, no bleeding, no downtime
- Compact tabletop unit · large touchscreen · continuous and pulse modes

TECHNICAL SPECIFICATIONS

Wavelength	980 nm
Power	30 W
Frequency	1–20 Hz
Spot	1.4 · 1.9 · 2.4 mm
Aiming beam	635 nm
Supply	100–240 V · 160 W



INCLUDED

Fibre handpiece, 2 tips and spot rings, installation, training and protocols

MANUFACTURER ADSS · EUDAMED-REGISTERED IMPORTER

WHY PHOIBOS TECHNOLOGIES

Certified devices

CE marked, with Diarmed Ltd registered as importer in the EU EUDAMED database.

Training at your premises

Hands-on operator training and treatment protocols for every device.

2-year warranty

Free servicing for two years, with maintenance contracts after.

Technicians island-wide

Technical support in Nicosia, Larnaca, Famagusta, Limassol and Paphos.

Flexible ownership

Purchase, instalments, rental or rent-to-own, to suit your business.

Free demo

Try any device at your premises, with no commitment.

WAYS TO ACQUIRE

Purchase

With a 2-year warranty and free service.

Instalments

Flexible repayment plan.

Rental

Monthly fee with service and maintenance.

Rent-to-own

Rent paid is deducted from the purchase price.

See it working, at your premises.

We bring the device, demonstrate it on a real client and show you the cost per session. No commitment.

+357 97 432 553

info@phoibostechnologies.com
phoibostechnologies.com

Diarmed Ltd · 169 Athalassas Ave, Office 303, 2015 Strovolos, Nicosia · Specifications may change without notice.